


Child Stool Collection Instructions (if your child is potty-trained)



1. Once your child is ready to poop, gather both child stool tubes from their packages, the spatula present in one package, and the toilet accessory.



2. Follow the instructions on the toilet accessory package to place the accessory on the toilet the child will be using before they sit down.



3. After your child has finished and you are ready, unscrew the purple tops from both tubes, while being careful to not spill the liquid inside the tubes.



4. Use the spatula to collect a pea-sized amount of stool that has not touched the toilet accessory or any other item.



5. Transfer the sample into the yellow tube top of one tube and repeat with another similar sized amount of stool for the other tube.



6. Use the spatula to level and make the sample flat, but do not push the sample into the tubes. Wipe the outside of the tubes with the provided cotton pad if needed.




7. Screw the purple cap back on tightly and shake the tube as hard as possible for 30 seconds. If large pieces of stool still remain, keep shaking.




8. Place the tubes into the *LabGuard* plastic bag and dispose of the other materials. Keep the samples in a room temperature location until they are ready to be shipped.

Child Stool Collection Instructions (if your child uses **diapers**)



1. Once your child has pooped in their diaper, set aside the diaper for stool collection.



2. Once you are ready, gather both child stool tubes from their packages and the spatula present in one package.



3. Unscrew the purple tops from both tubes, while being careful to not spill the liquid inside the tubes.



4. Use the spatula to collect a pea-sized amount of stool that has not touched the diapers, wipes, or any other item.




5. Transfer the sample into the yellow tube top of one tube and repeat with another similar sized amount of stool for the other tube. .



6. Use the spatula to level and make the sample flat, but do not push the sample into the tubes. Wipe the outside of the tubes with the provided cotton pad if needed.



7. Screw the purple cap back on tightly and shake the tube as hard as possible for 30 seconds. If large pieces of stool still remain, keep shaking.



8. Place the tubes into the *LabGuard* plastic bag and dispose of the other materials. Keep the samples in a room temperature location until they are ready to be shipped.